



EDGEAPEX HP™ EDGEFLOW™ EDGEPACK™ EDGEAPEX™

EDGEPOWER

MORE POWER FOR ENDODONTICS

The Edge Power portfolio was designed to work seamlessly with our **NiTi and Obturation solutions.**





EDGEAPEX HP™

EDGEAPEX HP™ Designed to enhance precision and efficiency in endodontic procedures. Cordless and with an Integrated Apex Locator for maximum convenience.

Benefits

- Compatible with the majority of the existing file systems for better adaptability.
- Exceptional versatility with two movement options, continuous rotation and reciprocation.
- Cordless portable motor with wireless charging, auto on/off, and an integrated Apex Locator for maximum convenience.
- High-performance brushless motor with quiet operation and an extended service life.
- It adopts real-time feedback technology and dynamic torque contra angle, effectively preventing file separation.

Specifications:

- Great flexibility with a 360° rotation contra angle.
- Ultra-small head and slim neck for enhanced visibility and improved access to difficult areas.
- 10 available memories and pre-program settings for EdgeEndo file systems to save time.
- 4 available motions: CW / CCW / REC / ATR
- Multiple Apical functions
- Right-handed and left-handed adjustable display
- RPM range: 100 - 2500
- Torque range: 0.4 - 5.0 Ncm



Item Code: **EAPEXHP**
 Motor Hand piece: **EAPEXRU**
 Charging Base: **EAPEXCB**

Accessories included in package

The box includes a motor handpiece, a charging base, and a contra angle. Measuring Wire, Lip Hook (x2), File Clip (x4), Touch Probe (x2), Power Adapter, Protective Silicon Cover (x2), Spray Nozzle, O-Rings (x2), Disposable Sleeves (1 pack).

Additional Accessories

Product	QTY per Pack	Item Code	Image
MEASURING WIRE	1 per Pack	EAPEXNMW	
TOUCH PROBE	2 per Pack	EAPEXTP	
LIP HOOK	2 per Pack	EAPEXLH	
CONTRA ANGLE HANDPIECE CA 161	1 per Pack	EAPEXCA	
FILE CLIP	4 per Pack	EAPEXFC	

Warranty: 2 years for handpiece motor, charging base, and power adapter. 1 year for contra angle attachment. 6 months for accessories.



EDGEFLOW™

360° Rotating with excellent gutta-percha fluidity for a uniform, and smooth, bubble-free output ensuring seamless fillings.

Benefits

- Excellent gutta-percha fluidity for a uniform, bubble-free output ensuring seamless fillings.
- Wireless and electric design for maximum comfort and ease of use.
- Dual battery system provides extended endurance, allowing for approximately 100 root canal fillings.
- Precise and intuitive display to better monitor the remaining gutta-percha.

Specifications:

- Auto heating stop after timeout; Auto shutdown after timeout
- Charging base with battery
- Temperature monitoring design
- Disposable needles come pre-packed with gutta-percha
- Reusable injection needles available in various sizes and lengths
- Preset temperatures: 100°C, 120°C, 150°C, 180°C, 200°C
- Injection Speeds: Slow, Medium, Fast



Accessories included in package

Item Code: **EEPTO**
Main Unit: **EEPTORU**
Charging Base: **EEPTOCH**

Main unit, Charging Stand, Power Adapter, Battery (x2), Heat Insulation Hood, Wrench of Needle, Gutta-Percha Ring, Injection Needle (x4) E23G-M (1) - E23G-L (1) - E25G-L (2).

Accessories included in package

Product	QTY per Pack	Item Code	Image
WRENCH	1 per Pack	EEPTOW	
LITHIUM-ION BATTERY	2 per Pack	EEPTOB	
THERMAL PROTECTOR CAP	2 per Pack	EEPTOHIH3	

Tips

Product	QTY per Pack	Item Code	Image
*E20G-S INJECTION NEEDLE	0 per Pack	EE20GSIN	
*E20G-M INJECTION NEEDLE	0 per Pack	EE20GMIN	
*E20G-L INJECTION NEEDLE	0 per Pack	EE20GLIN	
E23G-M INJECTION NEEDLE	1 per Pack	EE23GMIN	
E23G-L INJECTION NEEDLE	1 per Pack	EE23GLIN	
*E25G-M INJECTION NEEDLE	0 per Pack	EE25GMIN	
E25G-L INJECTION NEEDLE	2 per Pack	EE25GLIN	

Warranty: 1 year for device.

*only available upon request



EDGEPACK™

EDGEPACK™ Features high-precision temperature control with pre-set working temperatures, ensuring optimal compatibility with various gutta-percha types for precise and reliable obturation.

Features



- Wireless Gutta-Percha Obturation Device for seamless, cordless operation.
- Ergonomic and ambidextrous design for enhanced comfort and usability.
- 0.35mm pre-bend tip for improved access and better penetration into narrow root canals.
- High-precision temperature control for consistent and reliable obturation.
- Efficient GP cutting system with high-conversion system.
- Compact size for easy cleaning and convenient storage.
- Visual and auditory prompts for enhanced safety during procedures.

Specifications:


- Preset temperature options: 150°C, 180°C, 200°C, and 230°C.
- Advanced time-related features:
 - Ultra-fast warm-up time of 0.2 seconds.
 - Auto shut-off after 10 seconds for energy efficiency.
 - Long-lasting battery life—provides up to 4 hours per charge.

Accessories included in package









Item Code: **EEDPU**
Main Unit: **EEDDURU**

Charging Stand, Power Adapter, WP4004 Tip, WP4504 Tip, WP5506 Tip, WP5508 Tip, Heat Insulation Sleeve of Tips (x4).

Accessories included in package

Product	QTY per Pack	Item Code	Image
CHARGING BASE	1 per Pack	EEDPUCS	

Tips

Product	QTY per Pack	Item Code	Image	Product	QTY per Pack	Item Code	Image
WP4004 WORK TIP	1 per Pack	EEDPUT4004		EDGEPACK WP 5508 WORK TIP	1 per Pack	EEDPUT5508	
WP4504 WORK TIP	1 per Pack	EEDPUT4504		*EDGEPACK WP 5508L WORK TIP	0 per Pack	EEDPUT5508L	
*WP5004 WORK TIP	0 per Pack	EEDPUT5004		*EDGEPACK WP 5510 WORK TIP	0 per Pack	EEDPUT5510	
WP5506 WORK TIP	1 per Pack	EEDPUT5506		*EDGEPACK WP 6012 WORK TIP	0 per Pack	EEDPUT6012	

Warranty: 1 year for device and tips.

*only available upon request



EDGEAPEX™

EDGEAPEX™ Effectively locates the apical constriction, enabling precise determination of root canal working length for optimal treatment outcomes.



Item Code: **EEALU**
Main Unit: **EEALURU**
Charging Base: **EEALUB**

Accessories included in package

The EdgeApex is an apex locator and pulp vitality tester.



It utilizes cutting-edge multi-frequency technology combined with a highly efficient digital signal processor, delivering optimal accuracy and information processing speed.

Pulp vitality tests attempt to examine the presence of pulp blood flow to indicate a true health.

EDGEAPEX™ Apex Locator

- Charging Base (1)
- Measuring Wire (1)
- File Clip (4)
- Lip Hook (5)
- Touch Probe (2)
- Pulp Tester Probe (2)
- Tester (1)
- Lithium-ion Battery (1)
- Power Adapter (1)
- Endo Ruler (Available only upon request)
- Enlarged display of apical foramen (1)

Additional Accessories

Product	QTY per Pack	Item Code	Image
MEASURING WIRE	1 per Pack	EEALUMW	
TOUCH PROBE	2 per Pack	EEALUMF	
LIP HOOK	5 per Pack	EEALULH5	
ENDO RULER	Available only upon request	EEALUER	
FILE CLIP	4 per Pack	EEALUFC5	
PULP TESTER PROBE	2 per Pack	EEALUPTP	


Warranty: 1 year for device. 6 months for accessories.

EDGEPOWER WORKFLOW

Explore our complete range of endodontic products on our website: www.edgeendo.com

1. ACCESS

- EDGEAPEX HP™
- EDGE R GLIDE
- EDGE APEX™ Locator




2. SHAPING

- EDGEAPEX HP™
- EDGEOne R



3. IRRIGATION

- Paper Points



4. OBTURATION

- EDGEBIOCERAMIC™ Sealer
- Gutta Percha Points
- EDGECore™
- EDGEFill™
- EDGEFlow™
- EDGEPack™



FOLLOW US ON
SOCIAL MEDIA